



## Insulated Surgical Instruments

TSICS re-insulates Bi-polar and Mono-polar forceps, Vag. Speculums, Clip Applicers, L and J Hooks and many other instruments.

Process	Initial Insulation	Re-Insulation
Previous insulation stripping		●
Cleaning and degreasing	●	●
Primer added	●	●
Heating and powder coating applied in chosen colour	●	●
Electrical leakage testing	●	●

All instruments are cleaned, the insulation is re-tested for electrical leakage and proper operation, then is lubricated before being returned to the customer.

